


HAT Scanbody Kit

Einfach in der Anwendung | Höchste Präzision | Kosteneffizient
Das volle Potenzial von All-on-X-Behandlungen entfalten



Specification

Name	HAT Scanbody Kit
Accessories	Scanbody×12, Digital Lab Analog×24, Screwdriver×2
Compatibility	All Panda Scanner Series
Application	Edentulous Implantation
Material	Titanium Alloy, Medical-Grade Stainless Steel*
Surface Treatment	Sandblasted
Service Life	100 times

Supported Implant Systems

Supported Implant Systems

M1.4X0.3P

OSSTEM[®]
IMPLANT

straumann

Nobel
Biocare[™]

CORTEX[®]
The Future of Dental Implants

B&B DENTAL
IMPLANT COMPANY

ASTRA TECH
IMPLANT SYSTEM

DENTIS

MEGA[®]GEN
美格真

Dentsply
Sirona

Neo
Biotech
Dedication to Dentists

adin

HI[®]OSSEN
IMPLANT

AB

C-TECH
IMPLANT

BIOHORIZONS[®]
DENTISTRY • ORTHODONTICS • DENTURES

THOMMEN
Medical

CAMLOG[®]
SYSTEM

btk[®]
Implanting Trust,
in the right place

MICRO[®]CONE

NTA[®]IMPLANT
CONFIDENCE IN QUALITY

IEZ
IMPLANT Technology

NEODENT[®]
A Neodent Group Brand

S.I.N.
IMPLANT SYSTEM

IZEN
IMPLANT

J DENTAL CARE
just smile

Kontakt[™] PL
USER MANUAL

Anthogyr[®]
A Dentsply Sirona Brand

trau[®]gim 创英[®]

DentalMaster
百傲泰

ICX

SNUG
DENTAL IMPLANT

Bioconcept

M1.4X0.25P

BEGO

M1.6X0.2P

DIO[®] IMPLANT

M1.6X0.25P

antmed 安特医疗
DENTAL

Supported Implant Systems

M1.6X0.35P

Datsing 大清西格

WEGO 威高

MIS

Supported Implant Systems

M1.8X0.35P

Dentium 登腾

Cowellmedi
The Pioneers in Dental Implant and E-HBMF2

AlphaBio[™]
Implantology

