

ICX

The FAIR Premium Implant

NEW - ICX-MULTI® SUPERSTRUCTURES UP TO 65° ANGLE.

ALSO FOR ICX-ZYGOMA IMPLANTS.



**ICX-MULTI® WITH 65°
NOW AVAILABLE!**

Innovative and individual

determine and select the **ICX-MULTI®**.

Insert the ICX-MULTI® as usual in 0°, 17° or 35° and in gingival heights of 0.0mm, 0.5mm, 1.5mm and 3.5mm or now new and individual the ICX-MULTI® up to 65° angle and in the individual gingival heights 4.5mm up to 6.0mm

ICX-MULTI® – up to 65° Innovative and individual.



Use the freedom to create even complicated and difficult mastering prosthetic situations.

Increase your range of indications by using the new innovative and individual ICX-MULTI® superstructures.

These are compatible with all BoneLevel implants in the Diameters 3.75mm, 4.10mm and 4.80mm from medentis medical.

By synchronising the internal geometry, use all ICX-MULTI® abutments made of titanium with the ICX-PREMIUM, Use ICX-ACTIVE MASTER and ICX-ZYGOMA implants.

0° TO 65° ANGULATION



individual



ICX-MULTI® ABUTMENTS

Every new ICX-multi® design is individually developed for you for your needs created.

Order your individual ICX-MULTI® abutment today.
Vary up to 65° in individual Degrees and in the gingival heights 4.5mm up to 6.0mm.

Thus, from now on 53.5° GH 4.5mm or 63.0° GH 4,5mm possible. Choose freely and individually.

Important notice:

The angled, **individual** superstructures are with the 1.4 hexagonal instrument!



ALSO FOR ICX-ZYGOMA IMPLANTS

