

ICX

The FAIR Premium Implant



ICX-PREMIUM & ICX-ACTIVE MASTER TISSUE LEVEL

New!



Welcome to the ICX-Family!

Service-Tel.: +49 (0)2641 9110-0 · www.medentis.de

Mo.-Fr.: 7.30 bis 19 Uhr

medentis
medical

ICX-PREMIUM & ICX-ACTIVE MASTER TISSUE LEVEL

With the ICX-TISSUE LEVEL you now have the opportunity to react individually to the different bone conditions. Use the possibilities of the proven thread-designs with the advantages of a polished shoulder.

ICX-PREMIUM BL, ICX-ACTIVE MASTER BL as well as ICX-PREMIUM TL and ICX-ACTIVE MASTER TL have the same internal geometry. This means that you can continue to work with the system components work and take advantage of the modern ICX implant system.

Take advantage of the new short (6.5mm) and ultrashort (4.9mm) Ø 4.8 ICX implants.



ICX-BONELEVEL


ICX-TISSUE LEVEL

ADVANTAGES OF ICX-TISSUE LEVEL

- ✓ Fewer surgical interventions due to open healing
- ✓ One prosthetics line through platform switching, one surgery box
- ✓ Reduced plaque accumulation through machined surface
- ✓ Predictable aesthetics by shaping the soft tissue
- ✓ 1,7mm polished shoulder height (except for the Ø 4.8 x 4.9mm, here: 2.7mm tulip)

ICX-PREMIUM TL diameters and lengths:

Ø 3.45mm	in 6.5, 8, 10 and 12.5mm length
Ø 3.75mm	in 8, 10 and 12.5mm length
Ø 4.1mm	in 8, 10 and 12.5mm length
Ø 4.8mm	in 8, 10 and 12.5mm length



additionally

Ø 4.8mm	in 4.9mm and 6.5mm length!
---------	----------------------------


4.9mm length =
Deviating shoulder height

2.7mm



ICX-ACTIVE MASTER TL diameters and lengths:

Ø 3.75mm	in 8, 10 and 12.5mm length
Ø 4.1mm	in 8, 10 and 12.5mm length
Ø 4.8mm	in 8, 10 and 12.5mm length



ICX



**ICX-PREMIUM BL, ICX-ACTIVE MASTER BL
and
ICX-PREMIUM TL und ICX-ACTIVE MASTER TL
– only one prosthetic line!**

Welcome to the ICX-Family!

**Das FAIRE Premium
Implantat-System**

www.medentis.de

ICX

Fax to:

+49 2641 9110-120

Or by post: medentis medical GmbH · Walporzheimer Str. 48-52 · 53474 Bad Neuenahr/Ahrweiler

X	Article-No.		Quantity
		ICX-PREMIUM TL Implants	
	C-001-345065	ICX-plus ø 3.45 x 6.5mm ((Pack of 6)	
	C-001-345080	ICX-plus ø 3.45 x 8mm ((Pack of 6)	
	C-001-345100	ICX-plus ø 3.45 x 10mm ((Pack of 6)	
	C-001-345125	ICX-plus ø 3.45 x 12.5mm ((Pack of 6)	
	T-001-375080	ICX Premium TissueLine Implant 3.75 x 8mm ((Pack of 6)	
	T-001-375100	ICX Premium TissueLine Implant 3.75 x 10mm ((Pack of 6)	
	T-001-375125	ICX Premium TissueLine Implant 3.75 x 12.5mm ((Pack of 6)	
	T-001-425080	ICX Premium TissueLine Implant 4.1 x 8mm ((Pack of 6)	
	T-001-425100	ICX Premium TissueLine Implant 4.1 x 10mm ((Pack of 6)	
	T-001-425125	ICX Premium TissueLine Implant 4.1 x 12.5mm ((Pack of 6)	
	T-001-480049	ICX Premium TissueLine Implant 4.8 x 4.9mm ((Pack of 6)	
	T-001-480065	ICX Premium TissueLine Implant 4.8 x 6.5mm ((Pack of 6)	
	T-001-480080	ICX Premium TissueLine Implant 4.8 x 8mm ((Pack of 6)	
	T-001-480100	ICX Premium TissueLine Implant 4.8 x 10mm ((Pack of 6)	
	T-001-480125	ICX Premium TissueLine Implant 4.8 x 12.5mm ((Pack of 6)	
		ICX-ACTIVE MASTER TL Implants	
	TAM-001-375080	ICX Active Master TissueLine Implant Ø 3.75 x 8mm ((Pack of 6)	
	TAM-001-375100	ICX Active Master TissueLine Implant Ø 3.75 x 10mm ((Pack of 6)	
	TAM-001-375125	ICX Active Master TissueLine Implant Ø 3.75 x 12.5mm ((Pack of 6)	
	TAM-001-425080	ICX Active Master TissueLine Implant Ø 4.1 x 8mm ((Pack of 6)	
	TAM-001-425100	ICX Active Master TissueLine Implant Ø 4.1 x 10mm ((Pack of 6)	
	TAM-001-425125	ICX Active Master TissueLine Implant Ø 4.1 x 12.5mm ((Pack of 6)	
	TAM-001-480080	ICX Active Master TissueLine Implant Ø 4.8 x 8mm ((Pack of 6)	
	TAM-001-480100	ICX Active Master TissueLine Implant Ø 4.8 x 10mm ((Pack of 6)	
	TAM-001-480125	ICX Active Master TissueLine Implant Ø 4.8 x 12.5mm ((Pack of 6)	

Ku.-Nr./practice stamp/ surname/first name/postal code/place

medentis medical GmbH · Walporzheimer Str. 48-52 · 53474 Bad Neuenahr/Ahrweiler · Tel.: +49 (0)2641 9110-0

www.medentis.de